



Background

The BioFire FilmArray Pneumonia Panel *plus* identifies select organisms and antimicrobial resistance (AMR) genes in sputum (including endotracheal aspirate) and BAL (including mini-BAL) specimens using nested, multiplexed PCR followed by melting curve analysis.

AMR Gene Reporting

Viruses and atypical bacteria are reported qualitatively as Detected or Not Detected. The AMR genes are also reported qualitatively (Detected/Not Detected), but only if one or more applicable bacteria (i.e. potential carriers of the AMR gene) are also detected in the sample.

For bacteria assays, the BioFire Pneumonia Panel *plus* targets single copy genes and reports an approximate quantity of organism (i.e. bacterial DNA in copies/mL).

AMR Gene Result	Applicable Bacteria
<i>mecA/C</i> and <i>MREJ</i>	<i>Staphylococcus aureus</i> <i>Acinetobacter calcoaceticus-baumannii</i> complex
CTX-M	<i>Enterobacter aerogenes</i>
IMP	<i>Enterobacter cloacae</i> complex
KPC	<i>Escherichia coli</i>
NDM	<i>Klebsiella oxytoca</i>
OXA-48-like	<i>Klebsiella pneumoniae</i>
VIM	<i>Pseudomonas aeruginosa</i> <i>Proteus</i> spp. <i>Serratia marcescens</i>

Bacteria Reporting

Assay Result	Reported Result and Bin Result
Negative OR <10 ^{3.5} copies/mL	Not Detected
Positive AND ≥10 ^{3.5} – <10 ^{4.5} copies/mL	Detected 10 ⁴ copies/mL
Positive AND ≥10 ^{4.5} – <10 ^{5.5} copies/mL	Detected 10 ⁵ copies/mL
Positive AND ≥10 ^{5.5} – <10 ^{6.5} copies/mL	Detected 10 ⁶ copies/mL
Positive AND ≥10 ^{6.5} copies/mL	Detected ≥10 ⁷ copies/mL

Specimen Enrollment

A prospective clinical evaluation was conducted between October 2016 – July 2017 at 8 sites in the United States and enrolled a total of 1682 specimens (846 BAL and 836 sputum). Patient population included adults and children as well as hospitalized, out-patient, and ED.

Reference Methods

Pneumonia Panel	Reference / Comparator Method(s)
Bacteria	Method 1 – Sensitivity and Specificity performance: Quantitative reference culture (qRefCx) ^a
AMR Genes	Method 2 – Bin performance: quantitative PCR/sequencing method (qMol)
Atypical Bacteria and Viruses ^b	Conventional PCR & Sequencing

^a Specimens were plated at dilutions over four logs; bacteria considered “Detected” if enumerated at 10^{3.5} CFU/mL or greater

^b MERS-CoV assumed negative; not circulating in US during testing period

Performance: Viruses compared to Conventional PCR+seq

Analyte	Specimen Type	Panel Det.	Sensitivity/PPA TP/(TP+FN) %	Specificity/NPA TN/(TN+FP) %
Adenovirus	BAL	8	8/8 100	837/837 100
	Sputum	15	13/17 76.5	815/817 99.8
Coronavirus	BAL	31	18/21 85.7	810/823 98.4
	Sputum	34	28/32 87.5	796/802 99.3
Human Metapneumovirus	BAL	9	8/8 100	836/837 99.9
	Sputum	21	20/21 95.2	812/813 99.9
Human Rhinovirus/Enterovirus	BAL	63	52/54 96.3	771/782 98.6
	Sputum	109	96/96 100	717/730 98.2
Influenza A	BAL	13	10/10 100	830/833 99.6
	Sputum	16	13/13 100	819/822 99.6
Influenza B	BAL	6	5/6 83.3	837/838 99.9
	Sputum	14	12/12 100	821/823 99.8
Middle East respiratory syndrome coronavirus (MERS-CoV)	BAL	0	0/0 -	846/846 100
	Sputum	0	0/0 -	836/836 100
Parainfluenza Virus	BAL	18	16/18 88.9	824/826 99.8
	Sputum	30	28/29 96.6	804/806 99.8
Respiratory Syncytial Virus	BAL	3	3/3 100	841/841 100
	Sputum	47	43/43 100	787/791 99.5

Discrepancy investigation summary: Evidence of analyte presence was observed in 5/5 FN BAL, 10/10 FN sputum, 21/31 FP BAL, and 24/33 FP sputum specimens by various methods.

Performance: Bacteria compared to qRefCx

Analyte	Specimen Type	Panel Det.	Sensitivity/PPA TP/(TP+FN) %	Specificity/NPA TN/(TN+FP) %
<i>Acinetobacter calcoaceticus-baumannii</i> complex	BAL	7	0/0 -	839/846 99.2
	Sputum	28	10/11 90.9	807/825 97.8
<i>Enterobacter aerogenes</i>	BAL	13	6/7 85.7	832/839 99.2
	Sputum	12	3/4 75.0	823/832 98.9
<i>Enterobacter cloacae</i> complex	BAL	23	11/12 91.7	822/834 98.6
	Sputum	32	11/12 91.7	803/824 97.5
<i>Escherichia coli</i>	BAL	20	12/12 100	826/834 99.0
	Sputum	48	23/24 95.8	787/812 96.9
<i>Haemophilus influenzae</i>	BAL	82	10/10 100	764/836 91.4
	Sputum	107	16/18 88.9	727/818 88.9
<i>Klebsiella oxytoca</i>	BAL	11	2/2 100	835/844 98.9
	Sputum	19	9/9 100	817/827 98.8
<i>Klebsiella pneumoniae</i> group	BAL	27	15/15 100	819/831 98.6
	Sputum	65	21/23 91.3	769/813 94.6
<i>Moraxella catarrhalis</i>	BAL	29	0/0 -	817/846 96.6
	Sputum	75	5/5 100	761/831 91.6
<i>Proteus</i> spp.	BAL	9	5/5 100	837/841 99.5
	Sputum	23	15/15 100	813/821 99.0
<i>Pseudomonas aeruginosa</i>	BAL	74	36/36 100	772/810 95.3
	Sputum	160	103/106 97.2	673/730 92.2
<i>Serratia marcescens</i>	BAL	12	6/6 100	834/840 99.3
	Sputum	53	26/27 96.3	782/809 96.7
<i>Staphylococcus aureus</i>	BAL	116	46/47 97.9	729/799 91.2
	Sputum	204	111/112 99.1	631/724 87.2
<i>Streptococcus agalactiae</i>	BAL	25	1/1 100	821/845 97.2
	Sputum	43	9/9 100	793/827 95.9
<i>Streptococcus pneumoniae</i>	BAL	29	5/5 100	817/841 97.1
	Sputum	51	16/16 100	785/820 95.7
<i>Streptococcus pyogenes</i>	BAL	8	2/2 100	838/844 99.3
	Sputum	11	6/6 100	825/830 99.4

Discrepancy investigation summary: Evidence of analyte presence was observed in 2/3 FN BAL, 7/13 FN sputum, 326/328 FP BAL, and 545/547 FP sputum specimens by various methods. FN in 5 sputum specimens were due to misidentification of the organism by the comparator method, and 1 FN sputum was due to a specimen swap.

Performance: AMR genes compared to conventional PCR + seq (from specimen) when BioFire Pneumonia Panel *plus* had applicable detection

Analyte	Specimen Type	Panel Det.	Sensitivity/PPA TP/(TP+FN) %	Specificity/NPA TN/(TN+FP) %
CTX-M	BAL	6	6/7 85.7	144/144 100
	Sputum	9	8/10 80.0	280/281 99.6
IMP	BAL	0	0/0 -	151/151 100
	Sputum	0	0/0 -	291/291 100
KPC	BAL	3	2/2 100	148/149 99.3
	Sputum	7	7/7 100	284/284 100
<i>mecA/C</i> and <i>MREJ</i>	BAL	46	40/45 88.9	64/70 91.4
	Sputum	107	94/98 95.9	91/104 87.5
NDM	BAL	0	0/1 0	149/150 99.3
	Sputum	0	0/0 -	291/291 100
OXA-48-like	BAL	0	0/0 -	151/151 100
	Sputum	0	0/0 -	291/291 100
VIM	BAL	0	0/0 -	151/151 100
	Sputum	2	1/1 100	289/290 99.7

Discrepancy investigation summary: Evidence of AMR gene presence was observed in 7/7 FN BAL, 5/6 FN sputum, 6/8 FP BAL, and 13/15 FP sputum specimens by various methods.

Performance: Atypical Bacteria compared to conventional PCR + seq


Analyte	Specimen Type	Panel Det.	Sensitivity/PPA TP/(TP+FN) %	Specificity/NPA TN/(TN+FP) %
<i>Chlamydia pneumoniae</i>	BAL	1	0/0 -	844/845 99.9
	Sputum	0	0/0 -	835/835 100
<i>Legionella pneumophila</i>	BAL	2	2/2 100	833/833 100
	Sputum	0	0/1 0	826/826 100
<i>Mycoplasma pneumoniae</i>	BAL	4	3/3 100	841/842 99.9
	Sputum	7	7/8 87.5	827/827 100

Discrepancy investigation summary: *L. pneumophila* was observed by an independent molecular method in the single FN sputum specimen.

BioFire Pneumonia Panel *plus* Raw Value vs Quantitative Comparators


qMol was used to measure the quantity of bacteria present in a specimen. Individual specimen positions are determined by the qMol value on the y-axis and the BioFire Pneumonia Panel *plus* unbinned (raw) value on the x-axis. The region shaded in grey represents the range for values that span the 10^{3.5} Detected/Not Detected threshold and the highest value (10^{7.0}) validated by both the BioFire Pneumonia Panel *plus* and qMol comparator.

BioFire Pneumonia Panel *plus* unbinned values reported relative to the qRefCx are also presented to demonstrate how the BioFire FilmArray Pneumonia Panel result compares to a relevant gold standard.



BioFire Pneumonia Panel *plus* Bin vs Quantitative Comparators

The stacked area charts below demonstrate the proportion of each reported BioFire Pneumonia Panel *plus* bin result over the range of qMol and qRefCx values.



Bin Performance

BioFire Pneumonia Panel *plus* bin performance compared to qMol is shown for BAL and sputum. The qMol cells are broken into one-log ranges. A BioFire Pneumonia Panel *plus* bin result is considered concordant if the qMol value is within 0.5 log of the bin boundary. For example, the 10^{5.5} BioFire Pneumonia Panel *plus* bin (10^{4.5}-10^{6.5}) is concordant with the qMol range of 10⁴ to <10⁶. The top line (bold) in the table shows the qMol range broken down by one-log increments. The second line (in parentheses) shows the BioFire Pneumonia Panel *plus* bins that are concordant with that qMol range.

BAL Bin Performance						Sputum Bin Performance								
qMol Range (Concordant Bin)	ND to <10 ^{3.5} (ND)	10 ^{3.5} to <10 ^{4.0} (ND or 10 ^{3.5})	10 ^{4.0} to <10 ^{5.0} (10 ^{4.0} or 10 ^{4.5})	10 ^{5.0} to <10 ^{6.0} (10 ^{5.0} or 10 ^{5.5})	10 ^{6.0} to ≥10 ^{7.0} (10 ^{6.0} or 10 ^{6.5})	≥10 ^{7.0} (≥10 ⁷)	qMol Range (Concordant Bin)	ND to <10 ^{3.5} (ND)	10 ^{3.5} to <10 ^{4.0} (ND or 10 ^{3.5})	10 ^{4.0} to <10 ^{5.0} (10 ^{4.0} or 10 ^{4.5})	10 ^{5.0} to <10 ^{6.0} (10 ^{5.0} or 10 ^{5.5})	10 ^{6.0} to ≥10 ^{7.0} (10 ^{6.0} or 10 ^{6.5})	≥10 ^{7.0} (≥10 ⁷)	
Panel Bin	ND	12025	25	15	0	0	Panel Bin	ND	11392	51	47	6	1	
	10 ^{3.5}	47	28	36	5	0		10 ^{3.5}	79	31	84	18	2	
	10 ^{4.5}	5	6	47	45	5		1	10 ^{4.5}	12	9	75	80	20
	10 ^{5.5}	3	1	7	48	27		2	10 ^{5.5}	2	1	11	80	65
	≥10 ^{6.5}	2	0	0	17	74		66	≥10 ^{6.5}	4	0	0	13	113
% concordant		12025/12082 (99.5%)	53/60 (88.3%)	83/103 (79.0%)	93/115 (80.9%)	101/108 (93.5%)	66/70 (94.3%)		11392/11489 (99.2%)	82/92 (89.1%)	159/217 (73.3%)	160/197 (81.2%)	178/201 (88.6%)	207/224 (92.4%)
							786/931 (84.4%)							

Conclusions

- The BioFire Pneumonia Panel *plus* is a specific test for the detection of bacteria, atypical bacteria, viruses, and antimicrobial resistance genes from BAL and sputum specimens compared to molecular methods. Specificity of bacterial assays relative to culture was observed to be as low as 88% due to the fundamental differences between the molecular method used by the BioFire Pneumonia Panel *plus* (measuring target genomic copies) and culture (counting viable organisms), particularly for fastidious organisms.

- The BioFire Pneumonia Panel *plus* is a sensitive test for the detection of atypical bacteria, viruses, and antimicrobial resistance genes from BAL and sputum specimens. Sensitivity of bacterial assays relative to reference culture was excellent; five out of 16 total FN results in both specimen types were due to misidentifications made by the reference method, and one additional FN was attributed to a specimen swap.

- The binned bacterial result assigned by the BioFire Pneumonia Panel *plus* is highly concordant with an independent quantitative molecular method. However, the test consistently reports bacteria in sputum and BAL specimens at higher levels than reference culture. This is attributed to the fundamental differences between the molecular method used by the BioFire Pneumonia Panel *plus* (measuring target DNA copies) and culture (counting viable organisms), particularly for fastidious organisms.